

Supplier: Domel – European Union

CENTRIC CF108GR

MAIN ADVANTAGES:	<ul style="list-style-type: none">> High capacity> Easy loading (ergonomic design)> Thorough decontamination> Ease of use> End of run indication
TECHNICAL SPECIFICATIONS	
Model	CF108-GR Refrigerated Floor Standing
Power supply	230 V / 50 Hz 120 V / 60 Hz
Max. speed / RCF Swing-out rotor	4.500 rpm / 4.600 x g
Max. capacity Swing-out rotor	4 x 750 ml
Control system	microprocessor
Speed range	500 to 4.500 rpm
Speed precision	20 rpm
Temperature control range	-20°C to +60°C
Time	0 - 99h50 and continuous
Acceleration / Deceleration	10 / 10
Program settings	10 programs
Lid opening	motorized
Motor type	induction brushless motor
Certification	EC marked: IEC 1010-1; IEC 1010-2-020 directives
Dimensions (HxWxD)	893 x 536 x 671 mm
Crating (HxWxD)	1118 x 736 x 821 mm
Weight (Net/Crated)	130 / 145 kg