

CYTOLOGY SOLUTIONS

MAIN FEATURES

- > Adjusted for use in cytology, urology, microbiology, virology, oncology
- > Wide range of consumables for different applications
- > Centrifugation of up to 12 samples
- > Robust design
- > Microprocessor controlled functions of cover, time and speed
- > Adjustable acceleration and braking system
- > Special START - STOP features
- > FAST SPIN function
- > Possibility of programming

CENTRIC 350 Code 473.000	CENTRIC 250 Code 472.000	CYTO CENTRIFUGE
12	8	No. of places in rotor
2.500	2.500	Max. speed (RPM)
YES	YES	Brushless motor
49	48	Noise level (dBA)