

ROTAMIX SHP-10

STIRRER

MAIN FEATURES:

- > Small footprint design
- > Heating Features (SHP)
- > Heats to maximum temperature 300°C in 10 minutes

Stirring hot plate **ROTAMIX SHP-10** with heating plate $\phi 140$ mm is designed for stirring and heating liquids in laboratories.

ROTAMIX SHP-10 features a quiet motor drive with analog electronics.

Speed and temperature is fully adjustable.

Speed remains constant, even under varying loads.

Heating plate is protected with Teflon® coating.

- Ideal for stirring and heating applications
- Heats to maximum temperature in 10 minutes
- Max. capacity 10 liters
- Heating power 500 W
- Convenient, small footprint design
- Analog speed and temperature



TECHNICAL SPECIFICATONS

| | |
|-------------------------------|------------------------------|
| Model | RotaMix SHP-10 |
| Code no. | 690.000 |
| Voltage | 115V / 230V, 50 - 60Hz |
| Power | 540VA |
| Speed range | 200 to 1200 rpm |
| Stirring plate dimension | Teflon® $\Phi 140$, heating |
| Max. heating temperature | 300°C |
| Max. capacity | 10 liters |
| Ambient operating temperature | 4°C to 65°C |
| Net weight | 4,6kg |
| Dimensions | 170 x 265 x 135mm |

ORDERING INFORMATION

| | |
|----------------------------------|----------|
| ROTAMIX SHP-10, 230V / 50— 60 Hz | 744.400 |
| ROTAMIX SHP-10, 115V / 50— 60 Hz | 1744.400 |

We reserve the right to alter specification details etc. without prior notice or liability! V:11-19

DOMEL

DOMEL, d.o.o., Otoki 21, 4228 Železniki
BU Laboratory Systems, Na Plavžu 79, 4228 Železniki

Phone: +386 4 5117 525
E-mail: info.labs@domel.com
WEB: <http://labs.domel.com>

