



# MIXER B

REF : E051P-B

The **MIXER B** is an automatic balance able to shake, weight and control the flow of blood collection in real time.



- ✓ 100% traceability
- ✓ Blue LCD display
- ✓ Shows the collection status at real time by remote control
- ✓ Automatic and continuous check of the collection flow (ml/min)

- ✓ Automatic and manual tube clamping system for any type of tube with
- ✓ sensor to check the presence/absence of the tube in order to increase
- ✓ safety during blood collection
- ✓ Automatic and constant check of duration and speed of blood collection,
- ✓ collection flow and volume, date, time, alarms, etc.
- ✓ Tray suitable for any type of bags, removable and easy to clean
- ✓ Optimum security during blood collection
- ✓ Shaking speed easy to modify

# MIXER B

## TECHNICAL FEATURES

REF : E051P-B

<b>POWER</b>	Operating with either AC supply or rechargeable battery
<b>Main power:</b>	AC / DC
<b>Input:</b>	100 - 240 VAC / 47 - 63 Hz / 0,9 - 0,3 A
<b>Output:</b>	24 VDC - 1,25 A max
<b>Battery:</b>	Internal battery (up to 12 hours)
<b>STRUCTURE</b>	Thermoformed ABS
<b>Dimensions (cm) :</b>	31 x 29 x 17,5
<b>Packing (cm) :</b>	50 x 35 x 35
<b>Weight (kg) :</b>	3,5
<b>Gross weight (kg) :</b>	5,5
<b>Screen:</b>	Bleu LCD display easily visible from any angle
<b>Alarms:</b>	Visual and acoustic programmable alarms: <ul style="list-style-type: none"><li>- Duration and speed of the blood collection</li><li>- Codes control</li><li>- Low and high flow</li></ul>
<b>Shaking rate:</b>	Mixing cycle: 5 seconds (possibility to modify)
<b>Volume setting:</b>	Digital setting from 200 to 600 ml by intervals of 50 ml - Possibility to modify during donation - Automatic display conversion from grams to milliliters
<b>Accuracy level:</b>	+/- 1 ml
<b>Temperature:</b>	+5°C to 40°C
<b>Humidity:</b>	30% to 70%
<b>Storage conditions:</b>	-10°C to 50°C
<b>CERTIFICATIONS</b>	ISO 9001, ISO 13845, calibration Certificate

## USER RECOMMENDATIONS

**Calibration:** We recommend doing a tare (weighing) control at least once a month (see User's Manual)

**Cleaning:** Turn off the machine before cleaning. Clean with a damp cloth and do not use detergents.

## ACCESORIES



### Blood Bank Supplies

Via Faentina n. 207/Q – 48124 RAVENNA – ITALY T: +39  
0544 46404 F: +39 0544 501617

[www.vasinstrumenti.it](http://www.vasinstrumenti.it) [export@vasinstrumenti.it](mailto:export@vasinstrumenti.it)



**VASINI STRUMENTI S.r.l.**  
Divisione Elettromedicale