



## ARMARIO DE REFRIGERACIÓN LAB CARE PLUS / LAB CARE PLUS REFRIGERATOR



GAMA SERIES	Lab CARE PLUS	Lab CARE PLUS	Lab CARE PLUS	Lab CARE PLUS	Lab CARE PLUS	Lab CARE PLUS
MODELO MODEL	LTR17S	LTR17G	LTR25S	LTR25G	LTR30S	LTR30G
Rango T° Temperature range [°C]	+2 / +8	+2 / +8	+2 / +8	+2 / +8	+2 / +8	+2 / +8
Capacidad Nominal [litros] Nominal Capacity [litres]	170	170	250	250	300	300
Medidas Exterior [mm] L x P x A Exterior dimensions [mm] W x D x H	670 x 660 x 1285	670 x 660 x 1285	670 x 660 x 1530	670 x 660 x 1530	595 x 600 x 1868	595 x 600 x 1868
Medidas Interior [mm] L x P x A Interior dimensions [mm] W x D x H	570 x 523 x 800	570 x 523 x 800	570 x 523 x 1045	570 x 523 x 1045	488 x 463 x 1350	488 x 463 x 1350
Claro de puerta [mm] L x A Door opening [mm] W x H	568 x 610	568 x 610	568 x 855	568 x 855	483 x 1311	483 x 1311
Espesor Aislamiento [mm] Insulation thickness [mm]	50	50	50	50	55	55
Acabado Exterior Exterior finishing	Acero galvanizado plastificado blanco White plastic coated galvanised steel					
Acabado Interior Interior finishing	Acero galvanizado plastificado blanco White plastic coated galvanised steel					
N° Cajones Drawers Number	3	3	4	4	7	7
Dimensiones Cajón [mm] L x P Drawer Size [mm] W x D	500 x 445	500 x 445	500 x 445	500 x 445	420 x 400	420 x 400
Material Cajones Drawers Material	Acero galvanizado con pintura Epoxi Epoxy coated galvanised steel					
Intercambiabilidad de Estantes/Cajones Interchangeable Shelves/Drawers	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes
Tipo de Puerta Door Type	Ciega Solid	Cristal Glass	Ciega Solid	Cristal Glass	Ciega Solid	Cristal Glass
Cerradura Lock	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes
Salida USB USB Port	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes
Datalogger Datalogger	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes
Batería Battery	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional
Sonda simulación producto Product temperature simulation probe	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional
Sonda de temperatura ambiente Interior room temperature probe	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes
Termostato anticongelación Freezing protection controller	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional
Tipo de Controlador Controller type	Digital Plus	Digital Plus	Digital Plus	Digital Plus	Digital Plus	Digital Plus
Tipo de Refrigerante Refrigerant type	R134a	R134a	R134a	R134a	R134a	R134a
Potencia Eléctrica [W] Electric power [W]	271	280	305	315	420	436
Tensión/Frecuencia [V / Hz] Voltage/Frequency [V / Hz]	230 / 50	230 / 50	230 / 50	230 / 50	230 / 50	230 / 50
Peso [kg] Weight [kg]	91	96	110	118	119	124





LTR17





**170 L**

LTR25



**250 L**

LTR30



**300 L**


## ARMARIO DE REFRIGERACIÓN LAB CARE PLUS / LAB CARE PLUS REFRIGERATOR



GAMA SERIES	Lab CARE PLUS	Lab CARE PLUS	Lab CARE PLUS	Lab CARE PLUS	Lab CARE PLUS	Lab CARE PLUS
MODELO MODEL	LTR40S	LTR40G	LTR50S	LTR50G	LTR70S	LTR70G
Rango T° Temperature range [°C]	+2 / +8	+2 / +8	+2 / +8	+2 / +8	+2 / +8	+2 / +8
Capacidad Nominal [litros] Nominal Capacity [litres]	400	400	500	500	700	700
Medidas Exterior [mm] L x P x A Exterior dimensions [mm] W x D x H	670 x 660 x 1955	670 x 660 x 1955	687 x 705 x 2010	687x705x2010	687x815x2125	687x815x2125
Medidas Interior [mm] L x P x A Interior dimensions [mm] W x D x H	570 x 523 x 1470	570 x 523 x 1470	580 x 573 x 1542	580 x 573 x 1542	580 x 693 x 1529	580 x 693 x 1529
Claro de puerta [mm] L x A Door opening [mm] W x H	568x1280	568x1280	575x1348	575x1348	575x1488	575x1488
Espesor Aislamiento [mm] Insulation thickness [mm]	50	50	50	50	55	55
Acabado Exterior Exterior finishing	Acero galvanizado plastificado blanco White plastic coated galvanised steel					
Acabado Interior Interior finishing	Acero galvanizado plastificado blanco White plastic coated galvanised steel				Acero Inoxidable Stainless steel	
N° Cajones Drawers Number	8	8	8	8	8	8
Dimensiones Cajón [mm] L x P Drawer Size [mm] W x D	500 x 445	500 x 445	510x475	510x475	500 x 518	500 x 518
Material Cajones Drawers Material	Acero galvanizado plastificado blanco White plastic coated galvanised steel				Acero Inoxidable Stainless steel	
Intercambiabilidad de Estantes/Cajones Interchangeable Shelves/Drawers	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes
Tipo de Puerta Door Type	Ciega Solid	Cristal Glass	Ciega Solid	Cristal Glass	Ciega Solid	Cristal Glass
Cerradura Lock	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes
Salida USB USB Port	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes
Datalogger Datalogger	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes
Batería Battery	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional
Sonda simulación producto Product temperature simulation probe	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional
Sonda de temperatura ambiente Interior room temperature probe	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes	Si Yes
Termostato anticongelación Freezing protection controller	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional
Tipo de Controlador Controller type	Digital Plus	Digital Plus	Digital Plus	Digital Plus	Digital Plus	Digital Plus
Tipo de Refrigerante Refrigerant type	R134a	R134a	R134a	R134a	R134a	R134a
Potencia Eléctrica [W] Electric power [W]	305	325	651	671	432	456
Tensión/Frecuencia [V / Hz] Voltage/Frequency [V / Hz]	230 / 50	230 / 50	230 / 50	230 / 50	230 / 50	230 / 50
Peso [kg] Weight [kg]	133	141	150	158	153	160




LTR40




**400 L**

LTR50



**500 L**

LTR70



**700 L**


## ARMARIO DE REFRIGERACIÓN LAB CARE PLUS / LAB CARE PLUS REFRIGERATOR



GAMA SERIES	Lab CARE PLUS	Lab CARE PLUS	Lab CARE PLUS	Lab CARE PLUS
MODELO MODEL	LTR90S	LTR90G	LTR140S	LTR140G
Rango T° Temperature range [°C]	+2 / +8	+2 / +8	+2 / +8	+2 / +8
Capacidad Nominal [litros] Nominal Capacity [litres]	900	900	1400	1400
Medidas Exterior [mm] L x P x A Exterior dimensions [mm] W x D x H	757 x 978 x 2121	757 x 978 x 2121	1385 x 815 x 2125	1385 x 815 x 2125
Medidas Interior [mm] L x P x A Interior dimensions [mm] W x D x H	606 x 848 x 1512	606 x 848 x 1512	1277 x 693 x 1529	1277 x 693 x 1529
Claro de puerta [mm] L x A Door opening [mm] W x H	605 x 1494	605 x 1494	575 x 1488	575 x 1488
Espesor Aislamiento [mm] Insulation thickness [mm]	55	55	55	55
Acabado Exterior Exterior finishing	Acero galvanizado plastificado blanco White plastic coated galvanised steel			
Acabado Interior Interior finishing	Acero Inoxidable Stainless steel			
N° Estantes/Cajones Shelves/Drawers Number	8	8	16	16
Dimensiones Estante/Cajón [mm] L x P Shelf/Drawer Size [mm] W x D	570 x 790	570 x 790	500 x 518	500 x 518
Material Estantes/Cajones Shelves/Drawers Material	Parrilla de alambre de acero plastificado Plastic coated wire shelf		Acero Inoxidable Stainless steel	
Intercambiabilidad de Estantes/Cajones Interchangeable Shelves/Drawers	Si Yes	Si Yes	Si Yes	Si Yes
Tipo de Puerta Door Type	Ciega Solid	Cristal Glass	Ciega Solid	Cristal Glass
Cerradura Lock	Si Yes	Si Yes	Si Yes	Si Yes
Salida USB USB Port	Si Yes	Si Yes	Si Yes	Si Yes
Datalogger Datalogger	Si Yes	Si Yes	Si Yes	Si Yes
Batería Battery	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional
Sonda simulación producto Product temperature simulation probe	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional
Sonda de temperatura ambiente Interior room temperature probe	Si Yes	Si Yes	Si Yes	Si Yes
Termostato anticongelación Freezing protection controller	Opcional Optional	Opcional Optional	Opcional Optional	Opcional Optional
Tipo de Controlador Controller type	Digital Plus	Digital Plus	Digital Plus	Digital Plus
Tipo de Refrigerante Refrigerant type	R134a	R134a	R134a	R134a
Potencia Eléctrica [W] Electric power [W]	475	670	598	620
Tensión/Frecuencia [V / Hz] Voltage/Frequency [V / Hz]	230 / 50	230 / 50	230 / 50	230 / 50
Peso [kg] Weight [kg]	178	178	263	279




LTR90



900 L

LTR140



1400 L